

Bed Rails

Dropside Bed Rail



*Provides a support rail on the side of the bed
Ideal for patients requiring assistance while turning,
repositioning or during transfers
Easily installed and removed
Slides between mattress and bed base
Side rail drops to provide greater access
for bed transfers
Can be used for patient restraint
Powder coated mild steel construction*

Specifications:

Height: 900 mm
Width: 900 mm
Reach: 600 mm

HIRE: \$60 for 4 weeks

PURCHASE: \$240

Fixed Bed Rail



*Provides a support rail on the side of the bed
Idea for patients requiring assistance while
turning, repositioning or during transfers
Easily installed and removed
Slides between mattress and bed base
Powder coated mild steel construction*

Specifications:

Height: 480 mm
Width: 600 mm
Reach: 540 mm

HIRE: \$40 for 4 weeks

PURCHASE: \$65

Cobra Single Bed Pole



*Provides a support pole for the side of the bed
Ideal for patients requiring assistance while
turning, repositioning or during transfers
Slides between mattress and bed base
Can be positioned on either side of the bed
Easily installed and removed
Powder coated mild steel construction*

Specifications:

Height: 500 mm
Reach: 500 mm

PURCHASE: \$47.95

Self Help Pole Free Standing



*Free standing Self Help Pole assists user in
sitting or repositioning in bed
Sits behind or beside the bed
Handle attached to chain and is
height adjustable
Dismantles for storage or transportation
Powder coated mild steel frame construction*

Specifications:

Product Weight 7kg (base)
6kg (upright)
Max User Weight: 120kg

HIRE: \$40 for 4 weeks

PURCHASE: \$280

*NB: Terms and conditions apply to hired goods



GILLESPIES
HIRE & SALES SERVICE
"Enhancing Your Quality of Life"

Phone: 02 9411 2180
Fax: (02) 9411 2173
Email: sales@ghss.com.au
Web: www.ghss.com.au