

# Bedside Commodes



## STANDARD BEDSIDE COMMUNE

*Removable, padded seat with lift-up toilet seat and fixed bowl carrier beneath  
Raised Padded Armrests, removable bowl with lid and carry handle  
Removable back rest for easy storage or transportation  
Easy clean vinyl surfaces, three-sided shroud discreetly conceals bowl  
Powder coated steel frame that is strong, durable and easy to clean*

**Specifications:**

Overall width: 560mm  
Overall depth: 600mm  
Seat height range: 475-615mm  
Unit weight: 5.5kg  
Maximum User Weight: 124kg

**HIRE: \$60**  
for 4 weeks

**PURCHASE: \$220**



## BASKETWEAVE COMMUNE CHAIR

*A traditional wooden frame commode chair with woven fibre upholstery in white and beige  
The removable seat conceals a plastic commode pan and lid with carrying handle  
The high chair sides provide support to the user when rising or lowering*

**Specifications:**

Seat Height: 470mm  
Seat Width: 440mm  
Seat Depth: 440mm  
Backrest height: 775mm  
Weight: 7.2kg

**PURCHASE: \$630**



## DELUXE BEDSIDE COMMUNE

*Unique 'foldback' seat panel  
Padded backrest, seat and arms  
Shrouded/skirt panel  
Removable pan and lid*

**Specifications:**

Ajustable seat height 46-62cm  
Seat 46(w) x 50(d)cm  
Weight capacity 110kgs

**HIRE: \$80**  
for 4 weeks

**PURCHASE: \$399**

\*NB: Terms and conditions apply to hired goods



**GILLESPIES**  
**HIRE & SALES SERVICE**  
*"Enhancing Your Quality of Life"*

Phone: 02 9411 2180  
Fax: (02) 9411 2173  
Email: sales@ghss.com.au  
Web: www.ghss.com.au

13 Elizabeth Street, (PO Box 45), Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027  
Open Monday to Friday 9:00 am to 5:00 pm