

SEATING



Features:

*High back chair with lumbar and neck support
Adjustable legs
Adjustable seat depth
Weight capacity 120kgs*

Specifications:

Seat width 500mm
Seat depth 420-500mm
Backrest height 730mm
Adjustable seat height 425-575mm

HIRE: \$80 for 4 weeks

PURCHASE: \$525



Carequip 'Hunter' Chair

Features:

*Low back orthopaedic chair
Ideal for dining use or household activities
Adjustable seat height and angle provide optimum hip support
Curved padded backrest provides postural support and comfort
Comfortable padded armrests assists in standing and transfers
Powder coated steel frame
Maximum user weight: 100kg*

Specifications:

Seat: width: 440mm
Seat depth: 440mm
Seat height: 495-645mm
Arm height: 220mm
Back height: 450mm
Overall width: 580mm

HIRE: \$60 for 4 weeks

PURCHASE: \$299

*NB: Terms, conditions and deposits apply to hired goods

The 'Revolution' Chair

The revolution chair incorporates a unique, patented "swivel-slide" technology allowing easy access to the chair and to a table edge without moving its legs.

The lock/release mechanism allows the seat to gently slide back from the table so it can swivel left or right for easy exit.

Specifications:

Seat width: 450mm Overall width: 505mm
Seat height: 480mm Overall depth: 560mm
Seat depth: 485mm
Maximum user weight: 120kg

PURCHASE: Call for pricing



Made in Australia



Phone: 02 9411 2180
Fax: (02) 9411 2173
Email: sales@ghss.com.au
Web: www.ghss.com.au

13 Elizabeth Street, (PO Box 45), Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

Open Monday to Friday 9:00 am to 5:00 pm