

# HOSPITAL BEDS AND PATIENT LIFTERS

We hire and sell a range of bedroom equipment for both temporary & permanent requirements, such as remote controlled, electrically operated hospital beds and patient lifters to suit most needs. We also carry a range of pressure care overlays (please refer to our 'Pressure Care' brochure), and alternating air mattresses for those requiring additional relief.

Please feel welcome to call and discuss your specific needs with our experienced staff.

## K-DEE STANDARD BED



Four-section bed  
Easy to assemble  
Can be folded for easy transport  
Tilting  
Knee break

### Specifications:

Dual height adjustment (33-73cm or 40-80cm)  
Head and foot panels  
Four locking castors (4"/10cm)  
Linak Motors (Linak actuators)  
Maximum user weight: 180kg  
Adjustable back and foot rests (back 70°, foot 26°)  
Width: 440mm

## OXFORD ADVANCE FOLDING ELECTRIC HOIST

The Advance Hoist is a compact fold-away hoist yet it is able to perform a wide range of transfers. Whether from the floor, on to a bed, or wheelchair to commode chair, the Advance copes with ease. The Advance Hoist's unique 'swan neck' leg design allows the hoist to get even closer to wide obstacles. This is especially useful when encountering large chairs, bulky commodes and wheelchairs.

### Specifications:

Maximum working load: 155kg  
Maximum overall length: 1390mm  
Minimum overall length: 1290mm  
Width legs open external: 1360mm  
Overall height of legs: 110mm  
Front twin castors: 100mm  
Rear brake twin castors: 100mm

### Weights:

Mast, base and boom assembly: 28.9kg  
Power pack: 3kg  
Base assembly: 14.9kg  
Mast and boom assembly: 14kg

### Folded dimensions:

Length: 1180mm  
Height: 430mm  
Width: 550mm



**GILLESPIE'S** **TEL: 02 9411 2180**  
**HIRE & SALES SERVICE**

"Enhancing your quality of life"

Email: [sales@ghss.com.au](mailto:sales@ghss.com.au)  
Fax: (02) 9411 2173  
Web: [www.ghss.com.au](http://www.ghss.com.au)

13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027  
Open Monday to Friday 9am - 5.00pm

# Invacare Etude Plus Bed



## Features and benefits:

Built-in low brackets makes the Etude Plus height adjustable from 330mm-730mm or 400mm-800mm  
Choice of metal or wooden side rails available  
Transport kit available for easy storage and enhanced portability  
Range of accessories also include self-help pole, support handle and mattress support extension.

## Specifications:

Length outside: 2190mm  
Width outside: 1020mm  
Length inside: 2070mm  
Width inside: 900mm  
Mattress: 2000mm  
Net weight: 75.4kg  
Maximum User Weight: 180kg  
Mattress support dimensions:  
850-240-350-560mm

*Great value,  
Fast delivery!*



**GILLESPIE'S** TEL: **02 9411 2180**

**HIRE & SALES SERVICE**  
"Enhancing your quality of life"

Email: [sales@ghss.com.au](mailto:sales@ghss.com.au)

Fax: (02) 9411 2173

Web: [www.ghss.com.au](http://www.ghss.com.au)

13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027  
Open Monday to Friday 9am - 5.00pm