



**GILLESPIES**  
**HIRE & SALES SERVICE**  
*"Enhancing Your Quality of Life"*

# Pressure Care Catalogue

Quality patient pressure care products

Gillespie's stock a wide range of pressure care & patient comfort products to suit needs from low to high risk situations.

## Alleviating Constant Pressure

Please note, the information provided within is only a general guide to help with the prevention of sores, always seek professional medical advice.

- Keep patient off the sore area for periods of time
- Help the area breathe as much as possible
- Seek professional medical advice on how to best treat problem areas
- Reduce friction and sheen between patients skin and chair or bed

For more information on pressure care products or other equipment requirements, please contact our friendly staff or visit our website: [www.ghss.com.au](http://www.ghss.com.au)

Gillespies Hire & Sales Service  
13 Elizabeth Street, Artarmon  
Ph (02) 9411 2180  
Fax (02) 9411 2173  
[sales@gillespies.com.au](mailto:sales@gillespies.com.au)



**Foreword - Chris Gillespie**

Pressure reducing products help ease the threat of skin shear and skin friction on patients especially those with bony prominences.

Pressure ulcers are local loss of skin and soft tissue as a result of pressure from prolonged bed rest or sitting and are liable to form over any points at which pressure is exerted by the weight of the body. They can occur on the buttocks and on bony points such as heels, elbows, shoulders and the back of the head.

Anyone who sits or lies in one position for a long time might get much benefit from the products contained herein. The ability to widen the area of distribution for overall pressure and reduce the peak pressures are the aims of all the manufacturers of these products. Pressure sores can be serious. Consult your doctor if a sore is forming.



Kind Regards,



Chris Gillespie - Managing Director

**Store Information**

**PHONE:** 02 9411 2180  
**HOURS:** Monday to Friday 9.00am - 5.00pm  
**FACSIMILE:** 02 9411 2173  
**E-MAIL:** sales@ghss.com.au  
**WEB SITE:** www.ghss.com.au  
**ADDRESS:** 13 Elizabeth Street,  
 Artarmon NSW 2064  
**POSTAL:** PO Box 45, ARTARMON NSW 1570

**DELIVERIES:**  
 Our daily **metropolitan delivery service** is available to cater to those for whom visiting our store is either inconvenient or unnecessary. Delivery charges apply and are available on request.

*Contents*

*In Bed* ..... 03

*Pillows & Pressure Mattresses* ..... 04

*Seating Wedges* ..... 05

*Air, Gel Seat Cushions* ..... 06

*Lumbar Support* ..... 07

*Armchair Cushions* ..... 08

*Head, Heel & Elbow* ..... 09

*Medical Sheepskin* ..... 10

*Overlays & Medical Sheepskin* ..... 11

*Posture & Braces* ..... 12



**GILLESPIES**  
**HIRE & SALES SERVICE**  
*"Enhancing Your Quality of Life"*

**Contoured Bed Wedge - \$165**

Product Code - 3113



The Contoured Bed Wedge elevates your head and upper torso when sleeping or resting in bed to reduce or eliminate the symptoms of acid reflux or gastric reflux and other medical problems such as asthma, bronchial or breathing constraints, hernia and heart conditions.

Ideal for nursing homes, hospitals and home. Elevates and supports the head, neck, back and shoulders to support your body while sleeping or resting.

Pillow Slip White option \$32 or Waterproof \$45

**Leg Relaxer - \$140**

Product Code - LRQ



Designed for use on the bed to air circulation and vascular drainage  
Provides even weight distribution and alleviates lower back strain  
Ideal for those with varicose veins or swollen legs, phlebitis, pregnancy or convalesce  
Elevates the legs from the hips, taking the pressure off your lower back

**Pressure Care Assistant - \$106.30**

Product Code - SCPCA



The Shear Comfort<sup>®</sup> Pressure Care Assistant is a 65cm x 20cm strip of sheep skin with linen flaps which can be tucked in under the mattress to keep it in place. Areas at risk, such as the heels, sacrum, shoulders or head, can be individually cared for.

**Spine Reliever - \$249**

Product Code - 1671



A comfortable pillow for reclining in bed  
Provides the ultimate in comfort and support  
Designed for people with problem breathing, poor circulation, hiatus hernia, back or neck problems, acid indigestion or reflux  
Enables the user to get the body into a Zero-Gravity position, relieving stress and strain on the lower back

**BetterNeck Pillow - \$139**

Code - 1644



USA Made - insures the correct support for the neck and shoulder area. Relieves the cervical vertebrae, muscles and nerve paths. Use directly with supplied washable pillow case.

**Boomerang Pillow - \$60**

Code - 02BOOSTA



A comfortable tri-pillow for bed or chair. Contains 1700gms of blended mix fiber fill. Ideal for use when reading, etc. Quality filling provides comfort and support to the neck, back and shoulders.

**Overlay Superior - \$496.45**

Code - SCSO



A full covered bed overlay, designed to cover from head to toe. Suitable for long term patients. 3 - 9 full sheepskins are used to manufacture this quality product.  
**Size**  
Superior Overlay  
75cm x 150cm

**Overlay Silicon Fibre - \$380**

Code - 4523



Cylindrical hollow core fibres provide greater comfort and shear reduction. Soft breathable cotton top side. Water resistant polyurethane underside. Overlay can be reversed to suit. Elastic corner straps secure overlay in position. Machine washable.

**Ultimate 75™ - \$995**

Code - MK11



The Ultimate 75 is designed for home based patients or patients at moderate risk of pressure ulcers. End straps securely fasten the system to any single mattress.

**Ultimate 125 Mattress - \$1,725**

Code - MC-5-PA



Designed for patients of all risk levels up to a high risk of developing pressure ulcers. Featuring a Static button to activate full inflation of all cells enabling ease of process when transferring and treating patients, this premium 5 inch system is the ultimate in pressure ulcer treatment and prevention.

**Ultimate 200 Digital - \$2,450.00**

Product Code - PTM63540



The Ultimate 200 Digital replaces your existing mattress and is designed for patients of all risk levels up to a high risk of developing pressure ulcers. The whisper quiet digital unit provides ease of operation with a host of capabilities that ensure it is simple to operate and that the therapeutic performance is difficult to alter by accident. Featuring a Static button to activate full inflation of all cells enabling ease of process when transferring and treating patients, this premium 5 inch system is the ultimate in pressure ulcer treatment and prevention.

**Keyhole Cushion - \$82**

Code - KC



The Keyhole Cushion is a high density cushion with a symmetrical aperture. Designed for post natal or post operative discomfort, haemorrhoids, bruised coccyx and other painful conditions. The keyhole seat cushion gives therapeutic, pressure-free comfort when sitting.

**Coccyx Wedge - \$53.95**

Code - COX01



Promotes proper posture by enabling longer periods of sitting by supporting the back from underneath. It also has the added benefit of a cut out area relieving pressure on the tailbone and anal area.

**Seat Wedge Coccyx cut out \$59.95**

Code - 1641-BK



Helps improve spinal alignment for improved back posture. It is angled to neutralize the concavity of sunken and sloping seats that compress the body and spine.

**Backmate - \$189.95**

Code - BMEDBM



The go anywhere back support. It fully supports both the lower and middle regions of the back, reducing stressful loads on the spine.

**Events Cushion - \$19.95**

Code -CP4PH



Creates the ultimate comfort experience in the hardest of seats. Readily moulds to the shape of the body and seat, at the touch of a button. Easily deflated and folded for transportation.

**Diffuser Cushion - \$102**

Code - DIFF



A premium density pressure-relief cushion made from memory foam with a special inner core of heavy density traditional foam. Gently absorbs pressure and 'holds' the user securely and comfortably. Ideal as a wheelchair cushion or regular seating.

**Seat Wedge - \$69.95**

Code - 169



Designed to help relieve & prevent back pain, numb buttocks & other discomforts associated with long periods of sitting. Improves circulation in the legs. Posture in the car

**Cairecush - \$63.95**

Code - 3069



Lightweight PVC air cushion. Improves blood circulation to reduce risk of pressure related problems. Unique two layered design for greater immersion. Cell pattern improves surface stability. Easily cleaned with soap and water.

**Varilite Reflex Cushion - \$285**


Code - 73588



Varilite Reflex non-adjustable cushion is the choice for wheelchair users with a risk of tissue breakdown. Users and caregivers value the simplicity and ease of use of the Reflex. The ultra-lightweight design of the Reflex help conserve energy expenditure for clients who manually propel.

**ROHO Low Profile - \$752**

Code - IR109LPC



Cells are interconnected and their unique shape promote blood circulation to facilitate pressure prevention, healing and relief

Flexible cells adapt to patient contours and reduce shear forces while transferring or repositioning. As the body is immersed in the cushion, greater contact area is achieved for the dispersion of pressure, hence the pressure on any one area is minimised

**ROHO High Profile \$752**

Code - IR109C



Cells are interconnected and their unique shape promote blood circulation to facilitate pressure prevention, healing and relief Air pressure can be adjusted to suit individual requirements. Cushion comes fitted with a durable synthetic cover.

**Action Pilot - \$395**

Code - AK9008



Quality polymer gel which unlike fluid gel, does not displace from areas requiring pressure relief. Isolates and protects areas where shear forces and pressure points are relevant.

**Multi-Purpose Cushion - \$61**


Code - MP



Extremely durable dimplefoam which enables the bodyload to be gently spread. Ideal for anyone sitting for long periods. Available in corduroy or waterproof cover 42cm x 42cm x 5cm.

**Single Cell Archi - \$69.95**

Code - WCP1S



An innovative concept in seating and wheelchair cushions that addresses pressure distribution and positioning needs. Cushions mould to individual body contours by expelling air.

Also available in Dual Cell valve - \$269  
Code: WCP1D

**Ring Cushion 2 - \$79.95**

Code - RC



Made of natural Latex which because of its load bearing qualities needed for its donut shape is the ideal compound to provide adequate compression to withstand such pressures.

**Obusforme Seatpad - \$135**

Code - OBS



Lightweight polyurethane foam cushion. Assists shock and impact absorption. Ideal when sitting in a vehicle or when travelling. Encourages proper alignment of the pelvis and thighs, evenly distributing body weight for extended sitting comfort. Flexible sections to suit a wide range of seating surfaces.

**Obusforme Highback Backrest Support - \$152** Code -OBH **Lowback - \$135** Code -OBL



**Highback Features**

Added benefits and comfort of extended neck & head support. Ideal for those who spend long hours driving. Specially designed to accommodate vehicles with bench seating. Includes removable lumbar pad.

**Lowback Features**

Lightweight and portable, perfect for use at home, in the office or on the go! Includes a removable lumbar pad for additional back support that can be adjusted to fit your individual needs. Award winning design.

**Bifom Back Cushion - \$37** Code - 1663-BL




Promotes proper posture and helps alleviate discomfort associated with sitting for long periods. Includes zipped cover and straps to affix to chair.

**Lumbar Bucket Seat - \$55** Code - 1668-N



Ergonomically designed to promote proper posture for all day comfort, and back pain relief. Tapered back ideal for bucket seats in cars.

**Contoured Full back Support - Large \$65** Code - CBL **Med \$65** Code - CBM **Sm \$65** Code - CBS



Therapeutically designed to support the spine all the way! Carved from heavy duty foam, and tailored in beautiful hard wearing covers. Value price and ideal for car, couch and chair.

**Available in 3 sizes:**

- Small - 122cm to 158cm
  - Medium - 159cm to 186cm
  - Large - 187cm to 200cm
- 90% of people use medium**

**Lumbar Roll D - \$35** Code - LRD



Helps prevent the onset of lower back pain. Four sizes to hug all backs. D-shaped, Round, Oval and Egg. Helps maintain the spine's natural curvature.



D-Shape - Compact size, fits car seats and couches. Length 30 cm x 14 cm x 9 cm

**Lumbar support - \$15.95** Code - CP1PH



Provides support for lower back to help banish back pain. Ideal for use in a car, plane, train or bus. Expell air to mould to lower back Easily deflated for compact transportation.

**Back Eze I - \$59.95**

Code - 3092



**Back Eze II - \$71.95**

Code - 3093

The Original Moulded Back Support Cushion for relief, comfort and prevention of most lower back pain.

Improves posture, can improve health. Designed with a unique two way curve for optimum back support.

**Silicon Fibre 1 - \$115**

Code - CQ4521



Silicone coated fibre cushion, manufactured from cylindrical hollow core fibres for greater comfort and shear reduction. Soft breathable cotton top side, moisture resistant PVC underside. Machine washable.

**Silicon Fibre 2 - \$230**

Code - CQ4522



Silicone coated fibre filled wheelchair cushion, manufactured from cylindrical hollow core fibres for greater comfort and shear reduction. Back and arm pads for greater user comfort. Machine washable.

**Total Spine Support \$102**

Code - TS1



A full moulded back support featuring unique flanking support designed to prevent lateral movement. Relieves stress on the major structure of the lumbar spine and rib cage.

The perfect backrest support if you spend long hours driving.

**Wheelchair & Armchair cushions**

3 PRODUCTS



**Code - 90220 \$110**  
Wheelchair cushion with Techmaflex cover  
Size (cms) 9.5 x 40 x 40

**Code - 90221 \$120**  
Wheelchair cushion with Techmaflex cover  
Size (cms) 9.5 x 42 x 43

**Code - 90223 \$149.95**  
Armchair cushion with Techmaflex cover  
Size (cms) 9.5 x 50 x 50

The Techmaflex cover is made of a polyurethane coated polyester and features: Extendibility in all directions, soft steam-permeable material which regulates heat and humidity. Resistance to dust mites, bacteria and liquid. Anti static, washable, non slip underside, also fireproof.

**Head & Heel Holder - \$54** Code - 0735



Dual purpose pressure care item. As a **Heel Holder**, it is used to elevate the ankle away from the mattress, allowing air to circulate freely around the heel. This also allows a second person to change dressings easily.

As a **Head Holder**, it supports weak neck muscles whilst eating and reclining. It provides a sturdy support to ensure the neck is kept in a comfortable position.

**Heel & Elbow Protector - \$48** Code - 4524, 4525



Padded Heel or Elbow cushion. Provides comfort and assists in the prevention of pressure build-up in the heel or elbow. Internal seams are carefully positioned to avoid creating pressure points.

**Elbow Pads - \$45** Code - 4524  
One size



Provides comfort and assists in the prevention of pressure build-up in the elbow area. Manufactured from millions of silicone coated cylindrical hollow cored fibres. Silicone coating acts as a lubricant to eliminate friction and reduce shear forces. Adjustable velcro securing straps to suit individual requirements. Internal seams are carefully positioned to avoid creating pressure points.

**Prone Pillow - \$140** Code - 1921



Slightly raised with ventilation to the patient's face



Ventilation channel

This pillow will allow the patient to lie in the prone position with their head straight. The pillow is slightly raised with ventilation to the patient's face.

Commonly used in post eye surgery, to stop movement of the head.

The removable waterproof cover can be wiped clean or hand or machine washed, up to 80°C.

**Calf Plate Protector - \$47.05**

Code - SCCALF



Relieves pressure for wheelchair users.

Foot Plate Protector also available.  
Code - SCFOOT

**Elbow Protector - \$91.15**

Code - SCEPL



The Elbow Protector is designed to help the healing process as well as a preventative aid.

**Short Slipper Boot - \$149.95**

Code - SCSSBM  
Medium US Size 8-10



Ideal for indoor use and in bed where the foot is at risk of ulcers. The large tongue and generously sized Velcro straps allow for normal and swollen feet.

**Wrap Around Boot - \$134.75**

Code - SCWABM  
Size Medium 42-43



The toes are enclosed in this boot. With generous Velcro straps for easy closure and rubber soles for outdoor and indoor use.

**Tall Slipper Boot - \$POA**

Code - SCTSBL



Designed for vascular problems and is particularly suitable for those with diabetes and at risk of leg ulcers

**Overlay Superior - \$POA**

Product Code - SCSO



A full covered bed overlay  
Designed to cover from head to toe  
Suitable for long-term patients  
3-9 full sheepskins are used to manufacture this quality product  
Size: 75cm x 150cm

Please note:  
This is a special order item, lead time would be approximately 4-6 weeks.

### Convuluted Foam Mattress Overlay

Code - FNO24776  
Single Size



For low risk patients, soft open cell foam which provides additional comfort to existing mattress.  
Convuluted profile provides increased air circulation.  
Available in Single, Double or Queen sizes.

**Prices:**  
Single Overlay - \$151  
Double Overlay - \$170  
Queen Overlay - \$190

### Cushion it - \$124.75

Code - SCCI



Provides extra pressure care and comfort for standard seating and back products such as wheelchairs and standard foam seating.

### Sovereign Snug - \$105.85

Code - SBSSM



Beautifully crafted Snugs will accommodate foot oedema whilst offering support and comfort to the user. Made especially from carefully selected medical lamb-skin material.  
Sizes. XS, S, M, L



### Shear Comfort™ XD1900 Medical Overlay - \$154.05

Code - SCOL



Shear Comfort Medical Overlays offer superior pressure, shear, heat, and moisture dissipation for enhanced comfort and effective pressure ulcer prevention and treatment, helping avoid the discomfort, health risks, and financial burden of pressure ulcers and skin tears  
They can also be washed repeatedly at temperatures up to 80°C for effective infection control  
Shear Comfort XD 1900 is an advanced lambs wool product, boasting 1900gsm (typical woven wool products have 300gsm – 500gsm) for extra density, extra comfort, and extra care  
Can be used in the bed, when seated in a chair or used in the car  
Provides an interface between the patient hand surface.

### Heel Protectors - \$91.15

Code - SCHPM  
Size: Medium



Shear Comfort® Heel Protectors provide unencumbered security and protection for the at risk heel and ankle areas  
The thick wool redistributes pressure and reduces shear while the heat and water vapour dissipation qualities provide added skin protection and ensure comfort in use.

### Pressure Care Boot - \$134.75

Code - SCPCBS  
Size: Small 39-41



Shear Comfort® Pressure Care Boot features a rubber sole for use outside as well as in, plus an open toe for easy toe inspection, toe separation and air circulation.



**Shoulder Support - \$89.95**

Code - N71



Unique orthopaedic elastic support strap fits comfortably under the arm pit & assists to generate support to the shoulder. Assists soft tissue swelling. Self application & adjustment. Promotes greater blood circulation, Completely washable & hygienic..

**Hernia Support Belt - \$89.95**

Code - 590



Removable foam pads can be adjusted to suit left or right side. Made of lightweight material construction & assists to inhibit skin irritation. The belt can be worn whilst showering or swimming. In cases of post operative recovery it is recommended to remove both pads.

**Sacro Cynch Elastic Back - \$89.95**

Code - A495REG  
REG/M 67-105 cm



Made of lightweight lycra & orthopaedic elastic  
Delivers firm support  
Includes 2-layer closure  
First layer consists of lightweight elastic across the lower back, attached to lyra panels, delivering mild compression  
Outer straps deliver firm, powerful support when attached.  
The harder applied, the more compression  
Five pliable support stays

**Adjustable Elastic Ankle Brace - \$32**

Code - 3002  
Med. 45-50cm



This general purpose adjustable spiral lift elastic ankle brace, provides low level orthopaedic support & protection to painful, tired ankles  
Ideal as a support for sprained & strained muscles & ligaments.  
Includes an adjustable closure to assist in simple application & allows variation to compression levels  
Retains some body warmth & can increase blood circulation.

**Posture Shoulder Brace - \$49.95**

Code - 750



Assists with 'Round Shoulder' Syndrome. Easily applied, the brace improves posture position. The brace can be applied directly over bare skin with comfort. Improves posture, allows unrestricted movement, delivers all day comfort and is totally washable and hygienic.

**Deluxe Pouch Arm Sling - \$49.95**

Code - 861



Supports the arm during recuperation from broken bones, sprains and general injury. For maximum immobilization the Swathe secures the arm against the body. Fits left or right. Soft foam construction. Touch tape closure

**Deluxe Wrist Brace Thumb Spica - \$64.95**

Code - 282R2  
Medium Right



Palmar wrist splint & fully immobilized thumb splint in one  
For MP & CMC joints  
Indicated for Rheumatoid Arthritis & deQuervain's Syndrome  
Maximum contouring at the palmar crease, allows for full finger function  
Additional contouring at the fifth M.P.  
Allows for palmar spread  
Fully adjustable fit

**Elastic Hinged Knee Support- \$49.95**

Code - 43H2  
Medium 32-37



Rigid hinges minimize excessive medio-lateral movement of the knee joint, while offering excellent support and stability  
Made from lightweight medical grade orthopaedic elastic and features front closures for easier application  
Top & bottom straps lock the brace in position and minimise movement  
A perfect brace for many summer activities and when worn all day