

TOILET AIDS

Over Toilet Aid



Tubular Steel fabrication
Zinc treated frame to resist rusting
Adjustable Seat Height 480 - 610mm
Width - Arms 450mm / Base 540mm
Weight capacity 100kgs
Armrests to assist sitting & rising
Also in 'Wide/Heavy Duty' version

HIRE:
\$40 for 4 weeks

PURCHASE:
from: **\$115.00**

Toilet Surround



Tubular Steel fabrication
Phosphate treated to resist rusting
Adjustable Height 660 - 780mm
Width between armrests 450mm
Armrests to assist transfers

HIRE:
\$40 for 4 weeks

PURCHASE:
\$89.95

Raised Toilet Seat



Moulded plastic seat and lid
Available in 50mm, 100mm
& 150mm
Clamps to both sides
of toilet pan
Strong & lightweight
Weight capacity 190kgs
Simple to install and remove
Fits most toilet suites

HIRE:
from \$40 for 4 weeks
PURCHASE:
from **\$90**

Raised Toilet Seat with Arms



Designed to assist a person who has difficulty lowering themselves onto a standard height toilet seat. Designed to sit over existing bowl
Will fit all standard toilet seats
Foam Padded arms give exceptional support and safety and unit is easily transportable
Available in 50mm & 100mm

HIRE: \$40 for 4 weeks

PURCHASE:
\$115 50 mm
\$130 100mm

*NB: Terms, conditions and deposits apply to hired goods

GRAB RAILS:

We stock a range of steel grab rails which may be fitted to stable surfaces in the bathroom (walls & cabinets etc) to give support and aid stability. We supply specialised stainless steel rails, 'made to order' in various shapes to suit your particular requirements.

ACCESSORIES:

Also available are a wide variety of toilet & incontinence products including;

- Urinals (Male & Female)
- Bed Pans & Slipper Pans
- Mattress protectors (S,D,Q)
- Absorbant Bed Pads (Kylie / Drycare)
- Absorbant Chair Pads
- Blueys & Maxi Blueys

**Please feel free to call and discuss your specific needs with our friendly staff.
We look forward to being of service.**



GILLESPIES
HIRE & SALES SERVICE
"Enhancing Your Quality of Life"

TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au

13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

Open Monday to Friday 9am - 5.00pm