



GILLESPIES
HIRE & SALES SERVICE
"Enhancing Your Quality of Life"

Phone: 02 9411 2180
Fax: (02) 9411 2173
Email: sales@ghss.com.au
Web: www.ghss.com.au



Powerlift Recliner Ashley

Product Overview

Each chair in the Powerlift Recliner Ashley range is designed to suit a particular body shape and size. Each chair in the range is available in either one motor (back and legs elevate/recline simultaneously) or two motor (back and legs elevate/recline independently). The one motor version also features a unique wall saver function allowing the chair to recline whilst sitting in close proximity to a wall.

Quick Ship Colours available:

- ⊗ Top Suede: Mink (Brown), Lagoon (Green), Denim (Blue)
- ⊗ Cavatina: Solitude (Floral)
- ⊗ Lodge: Cypress (Velvet)

Powerlift Recliner Ashley Details

CODE	DESCRIPTION	RECOMMENDED USER HEIGHT	SEAT WIDTH	SEAT DEPTH	SEAT TO TOP OF BACK	WARRANTY	MAX. USER WEIGHT
CHP195900	Powerlift Recliner Ashley Petite 1 Motor	153 cm	440 mm	470 mm	580 mm	24 months	140 kg
CHP196000	Powerlift Recliner Ashley Petite 2 Motor	<153 cm	440 mm	470 mm	580 mm	24 months	140 kg
CHP195700	Powerlift Recliner Ashley Mini 1 Motor	153-160 cm	480 mm	470 mm	600 mm	24 months	140 kg
CHP195800	Powerlift Recliner Ashley Mini 2 Motor	153-160 cm	480 mm	470 mm	600 mm	24 months	140 kg
CHP195500	Powerlift Recliner Ashley Medium 1 Motor	163-178 cm	480 mm	520 mm	660 mm	24 months	140 kg
CHP195600	Powerlift Recliner Ashley Medium 2 Motor	163-178 cm	480 mm	520 mm	660 mm	24 months	140 kg
CHP195300	Powerlift Recliner Ashley Maxi 1 Motor	180-188 cm	530 mm	560 mm	720 mm	24 months	140 kg
CHP195400	Powerlift Recliner Ashley Maxi 2 Motor	180-188 cm	530 mm	560 mm	720 mm	24 months	140 kg
CHP196100	Powerlift Recliner Ashley Tall 1 Motor	191-198 cm	540 mm	570 mm	740 mm	24 months	140 kg
CHP196200	Powerlift Recliner Ashley Tall 2 Motor	191-198 cm	540 mm	570 mm	740 mm	24 months	140 kg