

# WALKING FRAMES



## Compact Seat Walker

Lightweight, compact and stable walker that is easy to manoeuvre that is easy to manoeuvre Easy to fold, making it transportable and features "Feather Touch" Handbrakes It also has a large seat with extra padding, as well as a shopping compartment under the seat makes this a versatile walker  
Floor to Seat Height: 480 mm  
Seat Width: 400 mm  
Maximum user weight: 100kg

**HIRE: \$40 for 4 weeks**  
**PURCHASE: \$235.00**



## Tri Wheel Walker

Power coated steel frame for long life Tri-wheel design with 8" wheels is easy to maneuver over any terrain Loop style brakes control speed and easily lock wheels  
Folds to only 9" for easy transport and storage  
Specifications:  
Width: 25.5"  
Length/Depth: 24.5"  
Height Range of Handles: 31.75"-36.75"  
Width Between Handles: 14.5"  
Width Between Rear Wheels: 21.5"  
Maximum user weight: 250lbs

**HIRE: \$40 for 4 weeks**  
**PURCHASE: \$199.95**



## Seat Walker with Hand Brakes

Seat Lightweight Walker  
Curved padded backrest  
Cushioned seat which lifts for access to the shopping basket  
Width - 620mm Depth 600mm  
Floor to seat height 550mm  
Seat width 440mm  
Safe working load 100kg  
Available colours:  
Red, Blue, Purple, Champagne

**HIRE: \$40 for 4 weeks**  
**PURCHASE: \$160.00**



## Lightweight Folding Aluminium Walker

Lightweight aluminium  
Height adjustable to suit individual requirements  
Folds easily for transport and storage  
Contoured plastic handgrips  
Fitted with anti-rattle bushes  
Maximum strength & stability  
Available in three different sizes  
Max user weight 127kg  
Front wheels and back skis available

**HIRE: \$40 for 4 weeks**  
**PURCHASE: \$110.00**  
**With Wheels/Skis: \$190.00**



## High Mack Bariatric Walker

Constructed from heavy-duty aluminium  
Hand Brakes, Indoor & Outdoor use  
Handle Height range 870 - 980mm  
Gross Weight 9.4kg - Capacity 150kg  
Base Size 700mm(w) x 640mm(d)  
Seat width 500mm  
Padded Seat & Removable Back Rest  
Front Swivel Castors, Folding Frame  
Front of Seat Carry Basket

**HIRE: \$100 for 4 weeks**  
**PURCHASE: \$360.00**



The ideal mobility solution for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot  
Front-wheel steering provides increased manoeuvrability  
Deluxe lever brake ensures safety  
Tool free height adjustment  
8" casters are ideal for indoor/outdoor use

**HIRE: \$100 for 4 weeks**  
**PURCHASE: \$450.00**

\*NB: Terms and conditions apply to hired goods



**GILLESPIES**  
**HIRE & SALES SERVICE**  
"Healthcare Equipment"

**TEL: 02 9411 2180**

Email: [sales@ghss.com.au](mailto:sales@ghss.com.au)

Fax: (02) 9411 2173

Web: [www.ghss.com.au](http://www.ghss.com.au)