

# WALKING FRAMES

## Compact Seat Walker



Lightweight, compact and stable walker that is easy to manoeuvre  
Easy to fold, making it transportable and features "Feather Touch" Handbrakes  
It also has a large seat with extra padding, plus a shopping compartment under the seat makes this a versatile walker  
Floor to Seat Height: 480 mm  
Seat Width: 400 mm  
Maximum user weight: 100kg

**HIRE: \$40 for 4 weeks**  
**PURCHASE: \$235.00**

## Aspire Vogue Indoor Walker



The stylish European designed Aspire Vogue Walkers are lightweight with a sturdy frame and easy folding system for compact transport.  
Smooth handle-height adjustment - Graduated and numbered shaft automatically locates to pre-set  
Mobile meal table – a user can then transport their meal using the tray on the stable walker to go to where they'd like to sit to have their meal.

Specifications:  
Length: 650 mm, Width: 540 mm  
Folded Width: 270 mm  
Bag SWL: 5 kg  
Tray SWL: 5 kg  
Max user weight 110kg

**HIRE: \$60 for 4 weeks**  
**PURCHASE: \$270.00**

## Seat Walker with Hand Brakes



Seat Lightweight Walker  
Curved padded backrest  
Cushioned seat which lifts for access to the shopping basket  
Width - 620mm Depth 600mm  
Floor to seat height 550mm  
Seat width 440mm  
Safe working load 100kg  
Available colours:  
Red, Blue, Purple, Champagne

**HIRE: \$40 for 4 weeks**  
**PURCHASE: \$160.00**

## Lightweight Folding Aluminium Walker



Lightweight aluminium  
Height adjustable to suit individual requirements  
Folds easily for transport and storage  
Contoured plastic handgrips  
Fitted with anti-rattle bushes  
Maximum strength & stability  
Available in three different sizes

Max user weight 127kg  
Front wheels and back skis available

**HIRE: \$40 for 4 weeks**  
**PURCHASE: \$110.00**  
**With Wheels/Skis: \$190.00**

## High Mack Bariatric Walker



Constructed from heavy-duty aluminium  
Hand Brakes, Indoor & Outdoor use  
Handle Height range 870 - 980mm  
Gross Weight 9.4kg - Capacity 150kg  
Base Size 700mm(w) x 640mm(d)  
Seat width 500mm  
Padded Seat & Removable Back Rest  
Front Swivel Castors, Folding Frame  
Front of Seat Carry Basket

**HIRE: \$100 for 4 weeks**  
**PURCHASE: \$360.00**



The ideal mobility solution for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot  
Front-wheel steering provides increased manoeuvrability  
Deluxe lever brake ensures safety  
Tool free height adjustment  
8" casters are ideal for indoor/outdoor use

**HIRE: \$100 for 4 weeks**  
**PURCHASE: \$450.00**

\*NB: Terms and conditions apply to hired goods



**GILLESPIES**  
**HIRE & SALES SERVICE**  
"Healthcare Equipment"

**TEL: 02 9411 2180**

Email: [sales@ghss.com.au](mailto:sales@ghss.com.au)

Fax: (02) 9411 2173

Web: [www.ghss.com.au](http://www.ghss.com.au)

13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

Open Monday to Friday 9:00am - 5:00pm