



GILLESPIE'S
HIRE & SALES SERVICE
"Enhancing your quality of life"

TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au

13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

Open Monday to Friday 9am - 5.00pm

Swift Wheelchair

Ultra
 Lightweight
 Only
 7.8 kg



Swift Wheelchair

The lightest folding wheelchair in its class

The Days Healthcare Swift has been designed to be a light, practical and easy to transport wheelchair, making it an ideal occasional use chair for both indoor and outdoor use.

Its compact dimensions allow good manoeuvrability indoors and the ultra light construction and ease to fold for storage or transport make it a very versatile chair.

Features:

- Ultra Light Wheelchair 7.8kg
- 115kg User Weight Capacity
- Folding Backrest (for storage or transport)
- Quick Release Wheels (for transport)
- Fixed Arms
- Swing away and removable leg rests
- Solid Tyres (no maintenance required)
- Extended brake handle (ease of brake application)

Specifications:

Weight:	7.8kg (frame only) 12.9kg (complete)
User Weight:	115kg
Seat Width:	460mm
Seat Depth:	420mm
Backrest Height:	410mm
Seat Height – Front:	500mm
Seat Height – Rear:	475mm
Armrest Height:	230mm
Front Castor Wheel:	8"x1" PU
Rear Wheel:	24" x 1 3/8" PU
Overall Width:	660mm
Upholstery :	Nylon (Flame Retardant)
Colour:	Silver



Maximum
 User Weight

Price: \$495

Hire: \$120*
for 4 weeks hire