




OMEGA TA1 WHEELCHAIR

Features:

- lightweight aluminium frame
- swing away removable arm rests
- swing away removable leg rests
- 12" rear wheels, attendant brake
- easy to fold for travel or storage
- user weight capacity of 125kg



OMEGA TA1 WHEELCHAIR	Metric	Imperial
Approximate Total Product / Unit Weight	11.1 kg	24.5 lbs
Overall Width (product or frame when opened)	570 mm	22.4 inch
Overall Depth (product or frame when opened)	930 mm	36.6 inch
Width / Depth (product or frame when closed)	280 mm	11.0 inch
Overall Height (or height range if adjustable)	920 mm	36.2 inch
Armrest Width (outside dimension)	560 mm	22.0 inch
Armrest Width (inside dimension)	450 mm	17.7 inch
Armrest Height (from top of seat)	210 mm	8.3 inch
Back / Back-Rest Height (from top of seat)	400 mm	15.7 inch
Seat Height (from top of seat to floor)	500 mm	19.7 inch
Seat Width (actual seat, not available space)	450 mm	17.7 inch
Seat Depth (actual seat, not available space)	410 mm	16.1 inch
Wheels (dimensions of front)	152 mm	6.0 inch
Wheels (dimensions of rear)	305 mm	12.0 inch
Recommended Maximum User Weight	125.0 kg	275.6 lbs

