



# GILLESPIE'S

## HIRE & SALES SERVICE

*"Enhancing your quality of life"*

Fact Sheet

# Sea Wheelchair

- Aluminium Frame
- Nylon Upholstery
- Retractable & Removable Armrests (height & depth)
- Removable Rear Wheel
- Height & Angle Adjustable Front Castor
- Retractable, Removable Footrests
- Adjustable Rear Axle for Optimal Positioning
- Adjustable Angle & Height Footplates
- Adjustable Tension Backrest & Seat
- Weight Capacity 120kg



Code: FL-0022a



**FONDLIGHT**  
AUSTRALIA  
Solutions for Life

Chair Height	Chair Length	Chair Width	Seat Depth	Seat Width	Seat Height	Rear Wheel	Front Castor	Chair Weight	Weight Capacity
900mm	1060mm	670mm	460mm	460mm	520mm	24 x 1 <sup>3</sup> / <sub>8</sub> "	8"	17.7kg	120kg

# 9411 2180

Monday - Friday 8:30 am - 5:00pm

Saturday - 9am - 12 Noon

13 Elizabeth St, Artarmon • [www.ghss.com.au](http://www.ghss.com.au)