


NITRO Rollators



The Nitro Rollator walking aid has been designed for easy manoeuvrability both indoors and outdoors. Featuring large 10" form casters, the Nitro has an enhanced turning radius thanks to the caster fork design.

NITRO ROLLATOR

- Lightweight aluminium frame which has been designed for easy manoeuvrability
- Suitable for use in tight spaces with a caster fork design that enhances turning radius
- Internally routed brake cables for added safety
- Push button height adjustment on handles to suit user requirements
- Cross brace design for added stability
- Easily folds with one hand to a compact size for storage and portability
- Integral seat handle serves as a useful handgrip for pulling rollator when used as a shopper or for carrying when folded
- Large 25cm (10") front casters for easy manoeuvrability
- Lightweight walking frame is suitable for indoor and outdoor use
- Removable zippered storage bag for storing shopping items and essentials
- Plastic clip holds the rollator closed when frame is folded for easy handling
- Back support height easily adjusts with tool free thumb screw
- Soft, flexible backrest provides comfort and security while seated
- Rear stepper tubes facilitate safe curb-climbing for both attendant and independent user



CODES: RTL10266AU - Red
RTL10266BLKAU - Black

Width	Depth	Height Min	Height Max	Seat Width	Seat Depth	Seat Height Min	Seat Height Max	Weight	Weight Capacity
58cm (23")	70cm (27.75")	71cm (28")	80cm (31")	46cm (18")	22cm (8.5")	46cm (18")	46cm (18")	8kg (17lb)	136kg (21st)

NITRO PETITE ROLLATOR

- Petite height ideal for users 5'2" and below
- Attractive, Euro-style design
- Brake cable inside frame for added safety
- Handle height easily adjusts with unique push button
- Back support height easily adjusts with tool-free thumb screw
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort
- Easily folds with one hand to ultra-compact size for storage
- Plastic clip holds the Nitro closed when frame is folded
- Cross-brace design allows for side-to-side folding and added stability
- Seat is durable and comfortable
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded
- Lightweight aluminium Frame

CODE: RTL10266-25AU - Red



Width	Depth	Height Min	Height Max	Seat Width	Seat Depth	Seat Height Min	Seat Height Max	Weight	Weight Capacity
59cm (23")	70cm (28")	71cm (27.5")	80cm (31")	46cm (18")	21.5cm (8.5")	46cm (18")	46cm (18")	8kg (17lb)	136kg (21st)