



OMEGA LA1 WHEELCHAIR

Features:

- lightweight aluminium frame
- swing away removable foot rests
- fold down back rest
- 9" rear wheels, attendant brakes
- easy to fold for travel or storage
- user weight capacity of 125kg



| OMEGA LA1 WHEELCHAIR | Metric | Imperial |
|--|----------|-----------|
| Approximate Total Product / Unit Weight | 8.9 kg | 19.6 lbs |
| Overall Width (product or frame when opened) | 640 mm | 25.2 inch |
| Overall Depth (product or frame when opened) | 1,060 mm | 41.7 inch |
| Width / Depth (product or frame when closed) | 290 mm | 11.4 inch |
| Overall Height (or height range if adjustable) | 940 mm | 37.0 inch |
| Armrest Width (outside dimension) | 560 mm | 22.0 inch |
| Armrest Width (inside dimension) | 440 mm | 17.3 inch |
| Armrest Height (from top of seat) | 220 mm | 8.7 inch |
| Back / Back-Rest Height (from top of seat) | 410 mm | 16.1 inch |
| Seat Height (from top of seat to floor) | 470 mm | 18.5 inch |
| Seat Width (actual seat, not available space) | 460 mm | 18.1 inch |
| Seat Depth (actual seat, not available space) | 400 mm | 15.7 inch |
| Wheels (dimensions of front) | 152 mm | 6.0 inch |
| Wheels (dimensions of rear) | 203 mm | 8.0 inch |
| Recommended Maximum User Weight | 125.0 kg | 275.6 lbs |

