



aspire

DA VINCI QUATTRO

aspire for... Mobility Support Comfort

Electric Lift & Recline Chair Features


- 4 x motor lift chair for recline, legrest, lumbar and headrest movements
- Hand control with integrated USB charging port
- Chaise padding design fully supports legs, redistributing pressure
- Available in a large range of colours to suit modern home environments
- Integrated headrest support allowing safe and accurate head positioning
- Integrated lumbar support promotes posture and seating position




FABRIC SWATCHES



Latte



Gunmetal



Fog



Dove

CODE	DESCRIPTION	SEAT WIDTH	SEAT DEPTH	SEAT-TO-FLOOR	SWL
CHP227500	Aspire Da Vinci - Quattro - Latte	508 mm	520 mm	495 mm	159 kg
CHP227510	Aspire Da Vinci - Quattro - Gunmetal	508 mm	520 mm	495 mm	159 kg
CHP227520	Aspire Da Vinci - Quattro - Fog	508 mm	520 mm	495 mm	159 kg
CHP227530	Aspire Da Vinci - Quattro - Dove	508 mm	520 mm	495 mm	159 kg

GILLESPIES

HIRE & SALES SERVICE

Phone: 02 94112180

Email: sales@ghss.com.au

Web: www.ghss.com.au

Unit 6, 3 Lanceley Place

Artarmon NSW 2064 AUSTRALIA

GHSS Pty Ltd - ABN 94 082 770 027