

# Bed Rails

## Dropside Bed Rail



*Provides a support rail on the side of the bed  
Ideal for patients requiring assistance while turning,  
repositioning or during transfers  
Easily installed and removed  
Slides between mattress and bed base  
Side rail drops to provide greater access  
for bed transfers  
Can be used for patient restraint  
Powder coated mild steel construction*

### Specifications:

Total Width: 900mm

**HIRE: \$80 for 4 weeks**

**PURCHASE: \$230**

## Fixed Bed Rail



*Provides a support rail on the side of the bed  
Idea for patients requiring assistance while  
turning, repositioning or during transfers  
Easily installed and removed  
Slides between mattress and bed base  
Powder coated mild steel construction*

### Specifications:

Height: 480 mm

Width: 600 mm

Reach: 540 mm

**HIRE: \$60 for 4 weeks**

**PURCHASE: \$65**

## Single Bed Pole, adj. crook handle



*Support pole designed to assist the user to  
balance while getting into and out of bed  
Adjustable crook handle can be located in  
four positions for added user comfort and security  
Base is positioned under the mattress  
Suitable for use on all bed sizes*

**HIRE: \$60 for 4 weeks**

**PURCHASE: \$85**

## Self Help Pole Free Standing



*Free standing Self Help Pole assists user in  
sitting or repositioning in bed  
Sits behind or beside the bed  
Handle attached to chain and is  
height adjustable  
Dismantles for storage or transportation  
Powder coated mild steel frame construction*

Max User Weight: 150kg

**PURCHASE: \$280**

**Please discuss your specific needs with our experienced staff**



**GILLESPIES**  
**HIRE & SALES SERVICE**  
*"Healthcare Equipment"*

**TEL: 02 9411 2180**

Email: [sales@ghss.com.au](mailto:sales@ghss.com.au)

Web: [ghss.com.au](http://ghss.com.au)

**Open Monday to Friday 9:00am to 5:00pm**

Aidacare Pty Ltd ABN 40 134 398 833

*\*NB: Terms and conditions apply to hired goods*

Unit 6/3 Lanceley Place Artarmon NSW 2064