

SEATING



Aspire High Back Classic Day Chair

A high back orthopedic chair that offers adjustability, comfort and support
 Contoured Backrest
 Softly padded armrests, seat and back upholstery are comfortable and supportive
 Easy clean vinyl surfaces
 Height Adjustable
 Sturdy Design – Generous 160kg weight capacity, yet lightweight
 Powder Coated Steel Frame
 Non-Slip Tips provide stability whilst sitting or transferring
 Available in either Slate or Champagne Vinyl

Specifications:

Seat width: 490mm
 Seat depth: 490mm
 Seat to floor height: 430-610mm
 Backrest height: 700mm
 Overall width: 630mm
 Unit weight: 16kg
 Maximum user weight: 160kg

HIRE: \$100 for 4 weeks

PURCHASE: \$450

Aspire Low Back Classic Day Chair



A versatile, low back orthopaedic chair
 Contoured Backrest
 Softly padded armrests, seat and back upholstery are comfortable and supportive
 Non-Slip Rubber Tips
 Easy clean vinyl surfaces
 Height Adjustable – Easy push button adjustment for comfortable positioning and changing needs
 Lightweight Design - lightweight yet sturdy design for easy manoeuvrability
 Powder Coated Steel Frame – Easy to clean, strong and durable
 Available in either Slate or Champagne Vinyl

Specifications:

Seat: width: 460mm
 Seat depth: 440mm
 Seat to floor height:: 430-605mm
 Back height: 449mm
 Maximum user weight: 160kg

HIRE: \$88 for 4 weeks

PURCHASE: \$359

Aspire Adjustable Day Chair



A high back orthopaedic chair that offers adjustability, comfort and support
 Depth Adjustable – Easy hand wheel adjustment for comfortable positioning and changing needs
 Contoured Backrest – High back design provides increased postural support
 Softly padded armrests, seat and back upholstery are comfortable and supportive
 Easy clean vinyl surfaces
 Height Adjustable – Easy push button adjustment
 Sturdy Design - Generous 175kg weight capacity
 Powder Coated Steel Frame – Easy to clean, strong and durable
 Non-Slip Tips provide stability whilst sitting or transferring
 Available in either Latte or Ink Vinyl

Specifications:

Seat Depth: 420-620mm
 Overall Depth: 760-860mm
 Seat Width: 510-610mm
 Backrest Height: 740mm
 Overall Width: 575-675mm
 Seat Height: 410-530mm
 Backrest Height: 740-790mm
 Overall Height: 1070-1220mm
 Unit Weight: 26kg
 Weight Capacity 160kg

HIRE: \$160 for 4 weeks

PURCHASE: \$809

Please discuss your specific needs with our experienced staff



GILLESPIES
HIRE & SALES SERVICE
"Healthcare Equipment"

Unit 6/3 Lanceley Place Artarmon NSW 2064

TEL: 02 9411 2180

Email: sales@ghss.com.au

Web: ghss.com.au

Open Monday to Friday 9:00am to 5:00pm

Aidacare Pty Ltd ABN 40 134 398 833

**NB: Terms and conditions apply to hired goods*